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Development of Drug Assay System Based on iPS Cell-Derived Neurons

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ABSTRACT

The effectiveness of therapeutic treatments for many neurological diseases remains extremely limited. Neurons generated from iPS cells are expected to contribute to the analysis of disease mechanisms and the development of therapeutic drug candidates. As part of Ricoh's biomedical business, we have developed technologies for drug evaluation using neuronal cells differentiated from iPS cells using unique differentiation techniques. Among these technologies, we describe recent progress in neuronal activity measurement with multi-electrodes, maturation of neurons, and analysis of neurons derived from Alzheimer's patients.

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
Fig.1 Co-culture of excitatory neurons and iPSC-derived astrocytes on MEA. Top: phase contrast image of neurons and astrocytes cultured in a well of a MEA plate. The 16 electrodes are visible. Bottom: Raster plot. The horizontal axis represents time. Each row shows firings recorded with an electrode. Green lines indicate synchronized firings. The continuous curve represents the rate of firing.

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Table 1 Compounds for drug response evaluation.

Compound	Receptor/ Ion channel	Mode	Concentration [μ M]	Neurons*
NMDA	NMDA receptors	Agonist	10, 100, 1000	EX
Glutamate	AMPA/ NMDA receptors	Agonist	0.3, 3, 30	GA
D-AP5	NMDA receptors	Antagonist	1, 10, 100	EX GA
GABA	GABA receptors	Agonist	0.1, 1, 10	EX GA
Picrotoxin	GABA receptors	Antagonist	0.1, 1, 10	EX GA
4-AP	Potassium Channel	Blocker	0.3, 3, 30	EX
CNQX	AMPA/ Kainate receptors	Antagonist	0.5, 5, 50	EX
Pentylenetetrazol (PTZ)	GABA receptors	Antagonist	10, 100, 1000	GA
Carbamazepine (CBZ)	Sodium Channel	Blocker	1, 10, 100	GA
Kainic acid	Kainate receptors	Agonist	0.1, 1, 10	EX GA
Dimethyl sulfoxide (DMSO)	-	Negative control	0.2%, 0.3%, 0.4%	EX GA

*EX: excitatory neurons, GA: GABAergic neurons


Fig. 2 Pharmacological responses of excitatory neurons co-cultured with iPSC-derived astrocytes. Each dot represents a well sample. Two-sided Wilcoxon rank-sum tests were performed to compare the effects of the drugs at each concentration with those at 0 μ M. *: $p < 0.01$.

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
Fig. 3 Pharmacological responses of GABAergic neurons co-culturing with iPSC-derived astrocytes. Each dot represents a well sample. Two-sided Wilcoxon rank-sum tests were performed to compare the effects of the drugs at each concentration with those at 0 μM . *: $p < 0.01$.

Excitatory neurons		
Drug	Conc. Range	Relative # of firings (\log_2)
NMDA	10, 100, 1000 μM	* *
D-AP5	1, 10, 100 μM	* * *
GABA	0.1, 1, 10 μM	* * *
Picrotoxin	0.1, 1, 10 μM	
Kainic acid	0.1, 1, 10 μM	
CNQX	0.5, 50, 500 μM	* *
4-AP	0.3, 3, 30 μM	* * *
DMSO	0.2, 0.3, 0.4%	

GABAergic neurons		
Drug	Conc. Range	Relative # of firings (\log_2)
Glutamate	0.3, 3, 30 μM	*
D-AP5	1, 10, 100 μM	*
GABA	0.1, 1, 10 μM	* * *
Picrotoxin	0.1, 1, 10 μM	* * *
PTZ	10, 100, 1000 μM	* * *
CBZ	1, 10, 100 μM	* *
Kainic acid	0.1, 1, 10 μM	* * *
DMSO	0.2, 0.3, 0.4%	

Fig. 4 Summary of pharmacological responses of Excitatory/GABAergic neurons co-cultured with iPSC-derived astrocytes. Firing number is shown as color-coded \log_2 fold change relative to the drug-free condition. Two-sided Wilcoxon rank-sum tests were performed to compare the effects of each drug with that of DMSO. *: $p < 0.01$.

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
Fig. 5 Advanced maturation of human iPSC-derived neurons.
 A: Schematic representation of the stages of neuronal development in iPSC-derived neurons generated with the Quick-Neuron™ technology (Ricoh Biosciences, Inc.).
 B: Time-course characterization of neuronal development in cultured iPSC-derived neurons with immunofluorescence staining of the axon-specific marker NF-H (green) and the dendritic marker MAP2 (red). Scale bar: 100 μm.
 iPSC: induced pluripotent stem cell; NF-H: neurofilament heavy chain; MAP2: microtubule-associated protein 2; and DAPI: 40,6-diamidino-2-phenylindole.

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
Fig. 6 Dendritic maturation and spine formation in iPSC-derived neurons.
 Drebrin clusters (green) accumulate in spine heads along the MAP2- positive dendrites (red).

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
Fig. 7 Correlations of the gene expression in iPSC-derived neurons with the developing human brain transcriptome (BrainSpan atlas⁹⁾). RNA-seq revealed the time-dependent transition from a prenatal (left, post-conception week 8-26) to a postnatal/adult-like profile (right, 4 months-40 years following birth).

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
Fig. 8 Differential expression of neuronal maturation markers in iPSC-derived neurons.

A: Expression of drebrin (all isoforms and isoform A)

B: Expression of tau (isoforms 3R and 4R)

Mean and standard deviation of the relative levels of expression determined by RT-qPCR from $n = 3$ biological replicates at each time point.

*: $p < 0.05$, ***: $p < 0.001$ from Dunnett's multiple comparison test for one-way ANOVA.

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
Fig. 9 AI-based imaging assay of the effect of glutamate on mature iPSC-derived neurons.

A: Non-treated control sample. White circles indicate examples of drebrin clusters identified by AI.

B: Cells treated with 100 μM glutamate for 10 min before immunostaining.

C: Dose-dependent response curve to glutamate.

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
Fig. 10 Spontaneous firing activity of excitatory iPSC-derived neurons measured by HD-MEA.

A: Firing rate distribution map across the electrode area over the indicated time course above. Scale bars: 500 μm .

B: Raster plots displaying detected spikes (dots) in 1,020 recording channels (y-axis) over 600 s (x-axis).

C: Network Assay results, showing the evolution of spontaneous network burst profiles over time.

D: Percentage of active electrodes across the culture area.

E: Number of network bursts per minute of recording.

F: Burst peak firing rate (average maximal firing rate during a network burst).

(n = 2 biological replicates)


Fig. 11 Dose-dependent response to synaptic blockers.
 A: AMPA/Kainate receptor antagonist NBQX.
 B: NMDA receptor antagonist D-AP5.


Fig. 12 Axon tracking assay.

- A: Maps of the spatial distribution of action potential amplitude and the reconstructed axonal propagation paths across the electrode area. Scale bars: 100 μm .
- B: Representative analysis of a single neuron.
- C: Action potential propagation after alignment of the signal traces along an axonal branch (red line in B).
- D: Linear regression interpolation to compute the conduction velocity.
- E: Average total axon length of single neurons.
- F: Average conduction velocity of axon branches.

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
Fig. 13 Pathology of Tau in Alzheimer's disease.

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The two iPS cell-derived excitatory neurons were Quick-Neuron™ Excitatory, Ricoh Biosciences, Inc.


Fig. 14 Accumulation of Tau in AD neurons.

iPSC-derived neurons at 1 to 8 weeks of culture were collected and the amount of Tau in the cell lysates was measured using the ELISA method. After 4 weeks of culture, the amount of Tau in AD neurons (CW50114) was higher than that of the healthy control (CW50065). $N = 7$, Error bar: SD. Student t-test ($*p < 0.001$).


Fig. 15 Accumulation of phosphorylated Tau in the cell body of AD neurons.

Neurons were stained for phosphorylated Tau (p-Tau, green) using fluorescence immunostaining. Blue Hoechst staining shows the nucleus of the cells. In AD (CW50114) neurons, p-Tau was accumulated in the cell body.

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- 1) Y. Shi et al.: Induced pluripotent stem cell technology: a decade of progress, *Nat Rev Drug Discov*, Vol. 16, pp. 115–130 (2017).
- 2) M. Yang et al.: Polyethyleneimine facilitates the growth and electrophysiological characterization of iPSC-derived motor neurons, *Sci Rep*, Vol. 14 (1), 26106 (2024).
- 3) Y. Ishibashi et al.: Verification of the seizure liability of compounds based on their in vitro functional activity in cultured rat cortical neurons and co-cultured human iPSC-derived neurons with astrocytes and in vivo extrapolation to cerebrospinal fluid concentration, *Toxicol Appl Pharmacol*, Vol. 476, 116675 (2023).
- 4) J. L. Wallace, A. A. Pollen: Human neuronal maturation comes of age: cellular mechanisms and species differences, *Nat. Rev. Neurosci*, Vol. 25, pp. 7–29 (2024).
- 5) Y. Ohara et al.: Early-stage development of human induced pluripotent stem cell-derived neurons, *Journal of Neuroscience Research*, Vol. 93, pp. 1804–1813 (2015).
- 6) G. Quadrato et al.: Cell diversity and network dynamics in photosensitive human brain organoids, *Nature*, Vol. 545, pp. 48–53 (2017).
- 7) A. Gordon et al.: Long-term maturation of human cortical organoids matches key early postnatal transitions, *Nature neuroscience*, Vol. 24, pp. 331–342 (2021).
- 8) W. Lin et al.: Dendritic spine formation and synapse maturation in transcription factor-induced human iPSC-derived neurons, *iScience*, Vol. 26, 106285 (2023).
- 9) M. Li et al.: Integrative functional genomic analysis of human brain development and neuropsychiatric risks, *Science*, Vol. 362, pp. 139–148 (2018).
- 10) K. Togo et al.: Postsynaptic structure formation of human iPSC cell-derived neurons takes longer than presynaptic formation during neural differentiation in vitro, *Molecular brain*, Vol. 14, 149 (2021).
- 11) L. Miguel et al.: Detection of all adult Tau isoforms in a 3D culture model of iPSC-derived neurons, *Stem Cell Research* 40, 101541 (2019).

- 12) L.S. Capano et al.: Recapitulation of endogenous 4R tau expression and formation of insoluble tau in directly reprogrammed human neurons, *Cell Stem Cell* 29, pp. 918–932.e8 (2022).
- 13) K. Hanamura et al.: High-content imaging analysis for detecting the loss of drebrin clusters along dendrites in cultured hippocampal neurons, *Journal of pharmacological and toxicological methods* 99, 106607 (2019).
- 14) K. M. LaBarbera et al.: Modeling the mature CNS: A predictive screening platform for neurodegenerative disease drug discovery, *Journal of neuroscience methods* 358, 109180 (2021).
- 15) T. Shirao et al.: The role of drebrin in neurons, *Journal of Neurochemistry*, Vol. 141, pp. 819–834 (2017).
- 16) Y. Sekino et al.: Activation of N-methyl-d-aspartate receptor induces a shift of drebrin distribution: Disappearance from dendritic spines and appearance in dendritic shafts, *Molecular and Cellular Neuroscience*, Vol. 31, pp. 493–504 (2006).
- 17) J. Müller et al.: High-resolution CMOS MEA platform to study neurons at subcellular, cellular, and network levels, *Lab on a Chip* 15, pp. 2767–2780 (2015).
- 18) S. Ronchi et al.: Hierlemann: Electrophysiological Phenotype Characterization of Human iPSC-Derived Neuronal Cell Lines by Means of High-Density Microelectrode Arrays, *Advanced Biology* 5, pp. 1–16 (2021).
- 19) B. Dejanovic, M. Sheng, J. E. Hanson: Targeting synapse function and loss for treatment of neurodegenerative diseases, *Nat Rev Drug Discov* 23, pp. 23–42 (2024).
- 20) O. D. Howes, E. C. Onwordi: The synaptic hypothesis of schizophrenia version III: a master mechanism, *Mol Psychiatry* 28, pp. 1843–1856 (2023).
- 21) H. Krasner et al.: From Stress to Synapse: The Neuronal Atrophy Pathway to Mood Dysregulation, *IJMS* 26, p. 3219 (2025).
- 22) C. F. Lippa et al.: A β -42 deposition precedes other changes in PS-1 Alzheimer's disease, *Lancet*, Vol. 352(9134), pp. 1117–1118 (1998).
- 23) J. Hardy, D. J. Selkoe: The amyloid hypothesis of Alzheimer's disease: progress and problems on the road to therapeutics, *Science*, Vol. 297(5580) pp. 353–356 (2002).
- 24) J. G. Wood et al.: Neurofibrillary tangles of Alzheimer disease share antigenic determinants with the axonal microtubule-associated protein tau (tau), *Proc Natl Acad Sci U S A*, Vol. 83(11), pp. 4040–4043 (1986).
- 25) M. T. Heneka et al.: Neuroinflammation in Alzheimer disease, *Nat Rev Immunol*, Vol. 25(5), pp. 321–352 (2025).
- 26) E. H. Corder et al.: Gene dose of apolipoprotein E type 4 allele and the risk of Alzheimer's disease in late onset families, *Science*, Vol. 261(5123), pp. 921–923 (1993).
- 27) J. Avila et al.: Role of tau protein in both physiological and pathological conditions, *Physiol Rev*, Vol. 84(2), pp. 361–384 (2004).
- 28) C. Despres et al.: Identification of the Tau phosphorylation pattern that drives its aggregation, *Proc Natl Acad Sci U S A*, Vol. 114(34), pp. 9080–9085 (2017).